

Multi-Drug Rapid Test Device (Urine)




Panel



Cup



Strip



Cassette

Intended Use

The Multi-Drug Rapid test Device (Urine) are rapid chromatographic immunoassays for the qualitative and simultaneous detection of the following drugs or drug metabolites in a variety of combinations in human urine. The designed cutoff concentrations and direct calibrator in urine at the following cut-off concentrations:

Test Item	Test Item	Calibrator	Cut-off (ng/mL)
ACE	Acetaminophen	Acetaminophen	5,000/1,000/500/300
ALC	Alcohol	Alcohol	20mg/dL
BAR	Barbiturates	Secobarbital	300/200
BUP	Buprenorphine	Buprenorphine	10/5
BZO	Benzodiazepines	Oxazepam	300/200/100
COC	Cocaine	Benzoyllecgonine	300/200/150/100
COT	Cotinine	Cotinine	200/100
DIA	Diazepam	Diazepam	300
EDDP	EDDP	2-Ethylidene-1,5-dimethyl-3,3-diphenylpyrrolidine	300/100
ETG	Ethyl Glucuronide	Ethyl Glucuronide	1,000/500/300
FYL	Fentanyl	Fentanyl	10/5
GAB	Gabapentin	Gabapentin	2,000
K2	Synthetic Marijuana	JWH-018, JWH-073 metabolite	100/50/20
K3	Synthetic	AB-Pinaca metabolite	20/10
KET	Ketamine	Ketamine	1,000/500
LSD	Lysergic	Lysergic acid diethylamide	50
MAM	6-mono-acetyl-morphine	6-MAM	10
MCAT	Methcathinone	Methcathinone	1,000
MDMA	Methylenedioxyamphetamine	d,l-Methylenedioxyamphetamine	500/100

Test Item	Test Item	Calibrator	Cut-off (ng/mL)
MDPV	3,4-methylenedioxyphenylpropanone	3,4-methylenedioxyphenylpropanone	1000
MEP	Mephedrone	Mephedrone	500
MET	Methamphetamine	d-Methamphetamine	1,000/500/300
MOP	Morphine	Morphine	300
MQL	Methaqualone	Methaqualone	300
MTD	Methadone	Methadone	300
NFYL	Norfentanyl	Norfentanyl	20/10
OPI	Opiate	Morphine	2,000
OXY	Oxycodone	Oxycodone	100
PCP	Phencyclidine	Phencyclidine	25
PGB	Pregabalin	Pregabalin	1,000/500
PPX	Propoxyphene	Propoxyphene	300
SOMA	Carisoprodol	Carisoprodol	1,000
TCA	Tricyclic Antidepressants	Nortriptyline	1,000
THC	Marijuana	11-nor-Δ9-THC-9 COOH	50/25
TRA	Tramadol	Tramadol	200/100
XYL	Xylazine	Xylazine	500
ZOL	Zolpidem	Zolpidem	50
ZOP	Zopiclone	Zopiclone	50
7-ACL	7-Aminoclonazepam	7-Aminoclonazepam	300

Product Highlights

- 1.Product Variety: The Multi-Drug One Step Screen Test Panel (Urine) offers any combination from 1 to 18 drugs of abuse tests.
- 2.High Quality : CE Approved.
- 3.Simple Procedure: Collect urine, dip, and read the results.
- 4.Fast Speed: Reading results at 5 minutes.
- 5.Highly Sensitivity and Specificity.
- 6.Convenient Storage: Room temperature storage.

operating steps

